



Product:

HAAKE Poly Lab QC
Twin Screw Extruder

Manufactured By:

Thermo Fisher Scientific

CTW100

Screw Diameter	Conical
L/D	-
Screw Set-up	3 Options
Gear Ratio	1 : 1
Rotating direction	Counter
Max. Screw speed	250 min ⁻¹
Max. Temperature	450°C
Max. Pressure	700 bar
Max. Torque	200 Nm
Heating Zones	3
Cooling	Air
External Heaters	2
Feeding Zone	Cooled
Additional feeding/Venting	-
Sensor Ports	2 x 1/2" UNF
Options	Backforce sensor, additional sensor ports, venting dome

Rheology Solutions is the sole Australian Distributor of this product. Full product information is available to download from www.rheologysolutions.com. If you have any questions or require further information please contact us. We look forward to being of further service.

Regards from the Team at Rheology Solutions.